InLabTech

Inhalation Equipment

INLABTECH JAPAN LLC 24-9-302 Wakamatsucho, Shinjuku-ku, Tokyo 1620056 Japan TEL. +81-3-4405-1701 FAX. +81-3-4405-1702
インラボテックジャパン合同会社 〒162-0056 東京都新宿区若松町24-9-302 TEL. 03-4405-1701 FAX. 03-4405-1702
e-mail: info@inlabtech.jp

<http://www.inlabtech.jp>

InLabTech

SAFE CONVENIENT OPERABILITY

We Strive to Bring New Possibilities in Your Research

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We design customized inhalation exposure systems for small animals such as, mice, rats, and guinea pigs. By choosing the number of cages and test substance generators, we can build a system that exactly meets the needs of your research. The small footprint of our system means there is no need to construct a special building dedicated to housing an inhalation exposure system. All our system requires is a corner in a laboratory or a small room with an appropriate ventilation duct.


* The products in this catalog are some of the recent custom made devices.

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INHALATION EXPOSURE SYSTEMS

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
MULTIPURPOSE INHALATION CHAMBER

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The MIC exposure chamber is an AAALAC standard compliant animal cage specially designed for the exposure purpose.

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Exposure Chamber Design

- The innovative top-to-bottom cyclone structure ensures the homogeneous distribution of aerosol throughout the chamber.
- The chamber is made of clear resin, and its high transparency makes for easy visual observation.
- Chemical-resistant plastic is used for the chamber material to assure a high level of reliability comparable to that of stainless steel.
- There can be no air leakage as the pressure inside the chamber is maintained slightly negative during the experiments.

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FEATURES

- The exposure chamber can hold small animals such as, mice, rats, guinea pigs, and rabbits.
- The number of exposure groups can be adjusted simply by adding or reducing the chambers.
- Daily cleaning is expedited due to the smaller chamber size.
- This whole-body exposure chamber can be used as a nose only exposure chamber with optional mouse holders.

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MIC UNIT

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The MIC Unit complies with OECD test guidelines and can be used to conduct any test levels by simply adding or removing units.

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N=5 (rat)
N=10 (mouse)
Dosage=6

FEATURES

- **Individual chamber control**
Each chamber is controlled individually, meaning that multiple substance-concentration setups are possible. This allows for simultaneous exposure experiments with varied test substance concentration levels.
- **Easy to use computerized recorder**
The MIC unit has a recording function, while its easy-to-read color display indicates real-time temperature, humidity, pressure, flow rate, and the concentration values of each chamber.
- **Various setup options**
By choosing the desired number of cages, test substance generators, and/or by combining MIC units, you can build a system that exactly meets the needs of your research. While the chamber by itself is for the whole-body exposure method, you can use mouse holders for nose-only exposure.
- **High-quality air treatment**
A dual system with an activated carbon filter and a HEPA filter ensures effective and efficient filtration of the contaminated air exhausted from the chambers.

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
The Compact Inhalation Exposure System is an all-in-one cabinet type apparatus, which consists of the computerized controller, air treatment system, and exposure unit. For the exposure unit, either Nose-only Holder Shelf Unit or an MIC can be chosen.

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FEATURES

- **A single unit for both whole-body and nose-only exposure methods**
Conventional inhalation exposure testing systems are difficult to configure to match changes in test plans and varying conditions because the choice has to be made at the time of installation between either a whole-body exposure system or a nose-only exposure system. With this system, the user can install a nose-only or a whole-body exposure chamber by themselves depending on the desired application so as to expand testing possibilities.
- **Ideal for experiments using a smaller amount of test substance**
The compact and well-organized air flow system enables efficient delivery of an aerosol generated from a smaller amount of test substance to chambers or holders. This means you can conduct experiments using smaller amount of valuable (and costly) test substances.
- **Air flow control**
Use of a digital mass flow controller ensures accurate flow control.
- **Data monitoring**
Aerosol concentration and data on environmental conditions (such as temperature, humidity, pressure, and flow rate) can be monitored and recorded in real time. The system offers accurate data that backs your experiments.
- **Equipped with an exhaust-treatment system**
Aerosol after exposure is treated using several devices including a cyclone dust collector, a bubbling filter, an activated carbon filter, and a HEPA filter.

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Application with Multipurpose Inhalation Cage (MIC)

The system can hold up to 5 rats or 10 mice in accordance with AAALAC standards. The aerosol distribution within the chamber is homogenized using a hydrodynamic design.

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The inhalation cylinder is compact and easy to handle.
The dual flow-past structure separates supply and exhaust air, and prevents gas exhaled by a test animal being inhaled by other animals.
As the total volume of the supply/exhaust air cylinder is small, a smaller amount of test substance is required.

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Touchscreen Controller
• Configures Supply and exhaust air flow parameters
• Monitor and record: temperature, relative humidity, chamber pressure, supply/exhaust air flow rate, aerosol concentration
Compliant with 21 CFR Part 11 (optional)

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Aerosol Distribution Tubes
Provide controlled amount of the aerosol to each chamber

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Integrated Pump and Mass Flow Controller
Provide precise control to the supply and exhaust air flow

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(From the left)
Bubbler, Cyclone Dust Collector, Pre-filter
Trap the aerosol efficiently after passing through the chamber

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Activated Carbon Filter and HEPA Filter
Treat the supply air before the exposure and the exhaust air after the exposure effectively

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Negative Pressure Isolator with two doors
• Prevents leakage of the test substance.
• Easy to install/remove the content through the front and the side doors

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
SIMPLIFIED INHALATION EXPERIMENT SYSTEM

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The Best Starter Package

Owing to its greater simplicity and operability, the system aids quicker implementation of inhalation exposure experiments. Its small footprint allows placement in an existing draft chamber, enabling you to achieve the desired experiment environment with minimum equipment investment.

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
This system was developed in collaboration with Division of Toxicology, Department of Pharmacology, Showa University.

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
The system can be readily combined with a wide variety of generators to expand the possibilities of your experiments.

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The latest data-processing program with precise aerosol control


Environmental monitoring items required as part of the experiment data are constantly monitored and recorded using precision sensors installed within the chamber. The system is in compliance with the FDA 21 CFR Part 11 standards required at GLP facilities.

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Aerosol after exposure undergoes disposal treatment before being discharged.

Aerosol after an exposure experiment passes through a bubbling cyclone dust collector (primary treatment) and then an activated carbon filter and a high-performance HEPA filter installed within the control panel before being discharged outside the room.

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The Dust Generator DG series is a line of vibration-free and noise-free solid microparticle generators perfect for the use in animal testing. This instrument is distinctive in its ability to generate small amount of microparticles at a consistent low concentration for a long period of time. The concentration accuracy is easily achieved by setting generation conditions on a touchscreen display.

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FEATURES

- **Consistent low-concentration generation**
Less than 1 gram of microparticles can be generated at a low concentration for an extended period of time.
- **Easy set up and data saving**
The generator conditions such as flow rate, metering ring velocity, and hopper feeder cycles are easily set up on a touchscreen. Up to five sets of conditions can be saved and loaded.
- **Automatic concentration control**
The combined use of a Dust Generator DG with a digital dust monitor enables automatic control of concentration inside an exposure chamber.
- **Real-time observation of the concentration**
By connecting a digital dust monitor, the real-time concentration changes can be observed on the generator touchscreen display. The data is saved on an SD card.

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The Mist Generator AG-5 can control precise amount of the test substance and the air flow rate to support consistent generation of the monodisperse aerosol. This single unit is suited to a wide range of applications from low flow rate in vitro tests to high flow rate in vivo tests.

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
The digital microsyringe pump sends the controlled and consistent amount of the test substance to the atomizer.

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The FlowMaster 50 with the integrated filter regulator gives precise control to the air pressure and flow rate.

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The specially designed two-fluid atomizer supports the consistent aerosol generation.

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The generation volume is adjustable by changing the internal orifice ring.

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Organic Solvent Gas Generator

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Organic solvent is heated and vaporized to generate gas for supply to each chamber with a concentration gradient. Generated gas can be supplied at a stable concentration.

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
The flow rate of the gas to the chamber is adjusted individually.

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The generator component is kept in a negative pressure to prevent the unexpected gas leakage.

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The generated gas is stabilized with a cooler.

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The organic solvent is bubbled and heated to generate gas.

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
Smoke Generator

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After igniting a combustible test substance using a stainless-steel pipe, the generated smoke can be supplied at a designated flow rate. This product is ideal for small-scale inhalation exposure experiments using combustible matter such as tobacco leaves and herbs.

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
VITROCELL® ALI Technology for Inhalation Toxicology

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
In response to the scientific need to expose in physiologically relevant conditions, VITROCELL® exposure modules have been specifically designed and engineered to enable direct exposure of mammalian cells or tissue at the air/liquid interface (ALI) where the cell systems are not covered with culture medium. Researchers can thus use all cell types cultivated on microporous membranes. This approach allows for more credible and authentic results than by submerged exposure due to a closer replication of the human physiology.

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VITROCELL® Exposure at The Air/Liquid Interface

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Phases	Phases	Phases
Gaseous	Gaseous	Gaseous
Semivolatle	Semivolatle	Semivolatle
Particle	Particle	Particle
A	A	A
B	B	B
C	C	C
Submerged Cultivation and Exposure in Incubator A Medium above cells B Cells on membrane C Medium below cells	Suspension Cultivation and Exposure in Incubator A Cells in medium	Air / Liquid Cultivation and Exposure in Exposure Module A Direct and controlled exposure of test atmosphere to cells B Cells on membrane C Medium below cells
インキュベーターでの液内培養とばく露 A. 細胞の上の培養液 B. メンブレン上の細胞 C. 細胞の下の培養液	インキュベーターでの浮遊培養とばく露 A. 培養液内の細胞	ばく露モジュールでのALI培養とばく露 A. 被験物質のコントロール及び細胞への直接ばく露 B. メンブレン上の細胞 C. 細胞の下の培養液
Interaction of test components with culture medium	Interaction of test components with culture medium	No losses No reaction of constituents with culture medium High sensitivity of system
Low sensitivity	Low sensitivity	
被験物質と培養液の接触がある → 培養液に被験物質が溶け込む 感度が低い	被験物質と培養液の接触がある → 培養液に被験物質が溶け込む 感度が低い	細胞のロスが無い 被験物質が培養液に溶け込まない 感度が高い

VITROCELL® Stainless Steel Cultivation and Exposure Module for Continuous Culture Medium Supply

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The VITROCELL® Stainless Steel module system is designed for the direct exposure of human cells to airborne substances such as gasses complex mixture, nano particles, and fibers. The system simulates the human exposure situation. The cultivation base module is built of stainless steel for maximum resistance.

Normally two modules are required: one module for exposure to the substances and a second module for clean air control.

The capacity can be easily increased by adding more modules. Each module can receive a different dose concentration, so that a complete dose-response profile is obtained in one experiment.

The cells are exposed at the air/liquid interface on 3 cell culture inserts using low flow rates of the aerosol.

After exposure, the cells are further processed to measure a wide range of endpoints, e.g. cytotoxicity, genotoxicity, proliferation, cellular and oxidative stress as well as inflammation.

Culture Medium Supply

The medium is supplied to the module by gravimetric method, syringe or using a medium pump for intermittent / continuous medium exchange.

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



The microbalance sensor can be fitted in the VITROCELL® Cloud, VITROCELL® 6 CF, VITROCELL® 12, VITROCELL® 12/6 CF, and VITROCELL® AMES Stainless Steel modules. It is capable of measuring the deposited mass in the module at a resolution of 10 nanogram/cm³ per second. After validation, they may be removed so that the experiments can be carried out using the cell culture insets. Alternatively, one microbalance sensor can remain in the module to monitor the dose during exposure. The remaining compartments of the module are used for cell culture inserts.

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With the VITROCELL® Microbalance Software, deposition data in ng/cm² is made attainable online. The data is logged into a .txt file, which is exportable as an Excel® file

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VITROCELL® Cloud

For the exposure to liquid aerosols

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VITROCELL® Cloud is specifically designed for dose-controlled and spatially uniform deposition of liquid aerosols on cells cultured at the air/liquid interface. This system can be used for aerosols generated from liquids and suspensions. The aerosol is applied for a short time of approx. 3–4 minutes.


The modules for 6 wells, 12 wells, and 24 wells are available. The aerosol particle size can be chosen from 4.0 μm – 6.0 μm or 2.5 μm – 6.0 μm . The cultivation module is made of electropolished stainless steel and kept at 37°C by a digital incubator. The Cloud aerosol chamber is made of Polycarbonate.

Possible fields of application include screening for inhaled drugs and toxicity testing of inhaled substances including nanoparticle suspensions.

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The smoking Robot VC 10® is specifically designed and manufactured to fulfill the requirements of whole smoke and aerosol generation for in vitro experiments. Suitable for conventional and electronic cigarettes, the VC 10® offers significant advantages over other commercial smoking machines:

Shortest distance from the smoke generation to the cell cultures
This system prevents smoke from aging and guarantee authentic smoke-composition when it reaches the exposure target.

Flexible system to incorporate other analytical tools
The robot is designed to allow easy access to all tubes, filters and the pumping system. Additional analytical equipment relevant to the experiment can be easily installed individually.


Freely programmable parameters
The computer system caters for highly flexible programming of the smoking process. All parameters of the smoking process such as puff duration, puff volume, puff frequency, and exhaust duration can be defined according to experiment requirements.

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In using a smoking machine, having to press the switch of a manual button activated e-cigarettes every 30 or 60 seconds is such a hassle. The VITROCELL Vapestarter is an automated solution which does the pressing work in a precise manner. The device can be connected to any VITROCELL smoking machines, and its trigger function is controlled by the software of the smoke generator.

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The system consists of Vapestarter Controller, e-cigarette holder and different Vapestarter units

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
Intratracheal instillation (ITI) of a toxic substance is an alternative administration method of inhalation exposure that requires complex aerosol generation, exposure chamber, Visual Instillobot Is an animal test device which was developed to perform the administration of IT/ more effectively, Also, Visual instillobot can be 100% success rate of administration.

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- Visual Instillobot
 - User friendly structure
 - Equipped with LED and CCD camera

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- My Com Controller
 - Adjustment of Injection amount

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- Foot switch
 - Stabilize instillation precision
 - Stable control of injection

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FEATURES


- High success rate
- Short instillation time (below 1 minute)
- No side effect of anesthesia and possible to quickly recovery
- Available repeated administration a day
- Administration of droplet or mist-spray

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
Application

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Evaluation of toxicity effect of PHMG-p by ITI : After ITI using Visual Instillobot, and then observations of biological anaylsis and histopathology

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Inhalation efficacy of Drug A in pulmonary fibrosis : After ITI using Visual Instillobot in pulmonary fibrosis animal model, observations of histopathology and BW change

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Option

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- Animal station is not lanrge, and there is convenient movement
- Possible to adjust the angle of holding plate of animal station, so possible to fit for each posture
- Also, possible to reduce the deviation of instillation depending on personnel difference
- Available using for mouse and rat

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Neurochemical Monitoring and Analytical Services

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FEATURES

Microdialysis and biosensor techniques allow in vivo monitoring of the release, uptake and metabolism of neurotransmitters and other neuroactive molecules diffusing in the extracellular space of the brain.

Studies including PK/PD and targeted bioavailability such as penetration of a drug through the blood-brain barrier could be carried-out in an elegant way by use of multiple microdialysis probes sampling from various body compartments.

Microdialysis offers the closest correlation of neurochemistry in vivo to behavior or to electrophysiological measures of brain activity.

Microdialysis strengthens the predictive validity of animal models of psychiatric and neurological/neurodegenerative diseases and serves as a unique tool for investigating the mechanisms of drug actions.

Microdialysis data may play a crucial role in the decision process concerning further development of a new drug.

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Pronexus Analytical AB is a contract research organization offering services in neurochemical in vivo monitoring and allied bioanalytical techniques. Our mission is to serve pharmaceutical companies and research organizations aiming to accelerate the process of drug discovery, evaluate mechanisms of drug actions and strengthen functional validation of candidate drugs. The bioanalytical services provided by Pronexus are focused on applications of techniques and methods for monitoring brain chemistry in vivo including microdialysis and biosensors. The bioanalytical laboratory and the animal research facility are GLP-compliant with full regulatory approvals.

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EXPERIMENTAL MODELS

Specialized skills in microsurgery, stereotaxic surgery and in vivo handling protocols

- Microdialysis protocols on anesthetized or awake rats and mice, proof-of-concept studies, monitoring neurotransmitters and other surrogate markers
- Pharmacokinetics and drug brain delivery, targeted drug delivery through the BBB, triple microdialysis – PK/PD profiles of unbound drug concentrations in blood and brain compartments
- Microsurgery, stereotaxic surgery on small rodents, injection/infusion cannulas including microdialysis probes, syringe minipump-based continuous/discontinuous drug infusions for up to 4 weeks, Alzet osmotic pumps
- Lesion models of neurodegenerative diseases.
- Simultaneous microdialysis sampling and behavioral monitoring

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ANALYTICAL SERVICES

Unsurpassed sensitivity of neurochemical HPLC methods

- Ultra-sensitive detection of monoamines
Unique fluorescence derivatization/HPLC methods developed for simultaneous determination of noradrenaline, dopamine and serotonin down to 10^{-11} M (100 attmol/20 µl) and for histamine (300 attmol/20 µl)
- HPLC techniques optimized for determination of neurotransmitters, neuro-modulators and related molecules in the microdialysates
Serotonin, noradrenaline, dopamine and their metabolites, acetylcholine, GABA, aspartate and glutamate, glycine, physiological amino acids, histamine, adenosine and other purines, cAMP, cGMP, kynurenic acid, radical oxygen species (ROS), $\text{NO}^{2-}/\text{NO}_3^-$, Mg^{2+} , Ca^{2+} and other ions
- HPLC and LC/MS methods for analysis of biological samples
Determination of neurotransmitters or specified biomolecules, drug levels in plasma, CSF and tissue extracts
- HPLC and enzymatic spectrophotometric assays for clinical microdialysates
Amino acid analysis, lactate, pyruvate, glycerol, glucose, enzymes and clinical biomarkers

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